

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\localizer

TA: 1:08 PAT: Off Voxel size: 2.2×1.1×10.0 mm Rel. SNR: 1.00 SIEMENS: gre

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 2	
Slices	20
Dist. factor	20 %
Position	Isocenter
Orientation	Transversal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Slice group 3	
Slices	1
Dist. factor	20 %
Position	Isocenter
Orientation	Coronal
Phase enc. dir.	R >> L
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	280 mm
FoV phase	100.0 %
Slice thickness	10.0 mm
TR	20.0 ms
TE	5.00 ms
Averages	1
Concatenations	26
Filter	Prescan Normalize, Elliptical filter
Coil elements	L4

Contrast

TD	0 ms
MTC	Off
Magn. preparation	None
Flip angle	40 deg
Fat suppr.	None
Water suppr.	None
SWI	Off
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	50 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved
Saturation mode	Standard
Special sat.	None
Tim CT mode	Off

System

Body	Off
L4	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1
Dark blood	Off
Resp. control	Off

Inline

Subtract	Off
Liver registration	Off
Std-Dev-Sag	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On
Wash - In	Off
Wash - Out	Off
TTP	Off
PEI	Off
MIP - time	Off

Sequence

Introduction	On
Dimension	2D
Phase stabilisation	Off
Asymmetric echo	Off
Contrasts	1
Bandwidth	180 Hz/Px
Flow comp.	No
Allowed delay	0 s
RF pulse type	Fast
Gradient mode	Normal
Excitation	Slice-sel.
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\AAHScout

TA: 0:28 PAT: Off Voxel size: 1.6×1.6×1.6 mm Rel. SNR: 1.00 SIEMENS: AALScout

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	Isocenter
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0 deg
AutoAlign	Head
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	260 mm
FoV phase	100.0 %
Slice thickness	1.6 mm
TR	3.15 ms
TE	1.37 ms
Averages	1
Concatenations	1
Filter	Normalize
Coil elements	BC

Contrast

Flip angle	8 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution

Base resolution	160
Phase resolution	100 %
Slice resolution	69 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	On
Intensity	Medium
Cut off	20
Width	4
Unfiltered images	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending

System

Body	On
LO1	Off
LO2	Off
LO3	Off
LO1	Off
LO2	Off
LO3	Off
LO1	Off
LO2	Off
LO3	Off

Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
Auto Coil Select	Off

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Inline

Time to center	11.5 s
----------------	--------

Sequence

Introduction	On
Dimension	3D
Asymmetric echo	Weak
Contrasts	1
Bandwidth	550 Hz/Px
RF pulse type	Fast
Gradient mode	Normal
Excitation	Non-sel.
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\me_mprage

TA: 8:07 PAT: Off Voxel size: 1.3×1.0×1.3 mm Rel. SNR: 1.00 USER: tfl_mgh_multiecho

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	R2.2 P1.7 F31.8
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	128
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	1.33 mm
TR	2530 ms
TE 1	3.31 ms
TE 2	6.99 ms
TE 3	8.85 ms
TE 4	10.71 ms
Averages	1
Concatenations	1
Filter	Normalize
Coil elements	BC

Contrast

Magn. preparation	Non-sel. IR
TI	1100 ms
Flip angle	7.0 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	75 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off

Normalize	On
Intensity	Medium
Cut off	20
Width	4
Unfiltered images	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Single shot
Series	Interleaved

System

Body	On
LO1	Off
LO2	Off
LO3	Off
LO1	Off
LO2	Off
LO3	Off
LO1	Off
LO2	Off
LO3	Off
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R2.2 P1.7 F31.8
Orientation	Sagittal
Rotation	0.00 deg
F >> H	256 mm
A >> P	256 mm
R >> L	171 mm

Physio

1st Signal/Mode	None
Dark blood	Off

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Off
Contrasts	4
Bandwidth 1	195 Hz/Px
Bandwidth 2	651 Hz/Px
Bandwidth 3	651 Hz/Px
Bandwidth 4	651 Hz/Px
Flow comp. 1	No
Flow comp. 2	No
Flow comp. 3	No
Flow comp. 4	No
Echo spacing	13.8 ms
RF pulse type	Fast
Gradient mode	Fast
Excitation	Non-sel.
RF spoiling	On
Readout polarity	Positive
Readout trajectory	Bipolar
Add. scale factor	8.0
Gradient spoiling	Integral
Gradient moment factor	3.0
Siemens reconstruction	On
Save raw k-space data	Off
Averaging	RMS

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\epi_slicepositioning

TA: 2.2 s PAT: Off Voxel size: 3.5×3.5×3.5 mm Rel. SNR: 1.00 SIEMENS: ep2d_pace

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	37
Dist. factor	10 %
Position	L0.0 P11.2 H13.8
Orientation	Transversal
Phase enc. dir.	P >> A
Rotation	-180.00 deg
Phase oversampling	0 %
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.5 mm
TR	2180 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	75 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	64
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Prescan Normalize	Off
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On

Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	Head > Brain
Auto Coil Select	Default

Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 P11.2 H13.8
Orientation	Transversal
Rotation	-180.00 deg
R >> L	224 mm
A >> P	224 mm
F >> H	143 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active
Meas[19]	Active
Meas[20]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	3004 Hz/Px
Free echo spacing	Off
Echo spacing	0.43 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres

TA: 3:05 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	19
Dist. factor	0 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	151 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	1500 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	120
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	107 mm
Position	L0.0 P32.3 F24.8
Orientation	C > T37.4
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Rotation	180.00 deg
R >> L	151 mm
A >> P	76 mm
F >> H	35 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres_shorttr_stack

TA: 3:15 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	0 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	400 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	480
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	107 mm
Position	L0.0 P32.3 F24.8
Orientation	C > T37.4
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	9 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres_shorttr_stack

TA: 3:15 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	0 %
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	400 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	480
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	110 mm
Position	L0.0 P22.0 F18.1
Orientation	C > T39.5
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	9 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres_shorttr_stack

TA: 3:15 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	0 %
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	400 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	480
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	110 mm
Position	L0.0 P22.0 F18.1
Orientation	C > T39.5
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	9 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres_shorttr_stack

TA: 3:15 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	0 %
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	400 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	480
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	110 mm
Position	L0.0 P22.0 F18.1
Orientation	C > T39.5
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	9 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_highres_shorttr_stack

TA: 3:15 PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 SIEMENS: ep2d_bold

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	On
Start measurements	single

Routine

Slice group 1	
Slices	5
Dist. factor	0 %
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	400 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	70 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	480
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	7/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	110 mm
Position	L0.0 P22.0 F18.1
Orientation	C > T39.5
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A58.9 H21.3
Orientation	T > C-39.0
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	9 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	15
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

SIEMENS MAGNETOM TrioTim syngo MR B17

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\gre_field_mapping

TA: 4:20

Voxel size: 2.3×2.3×5.0 mm

Rel. SNR: 1.00

SIEMENS: gre_field_mapping

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
FoV read	300 mm
FoV phase	100.0 %
Slice thickness	5.0 mm
TR	1000 ms
TE 1	10.00 ms
TE 2	12.46 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BC

Contrast

MTC	Off
Flip angle	90 deg
Fat suppr.	None
Averaging mode	Short term
Reconstruction	Phase
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Special sat.	None

System

Body	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Rotation	90.00 deg
A >> P	300 mm
R >> L	300 mm
F >> H	5 mm

Sequence

Introduction	On
Dimension	2D
Asymmetric echo	Off
Contrasts	2
Bandwidth	260 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\vessels_head

TA: 0:15

Voxel size: 1.5×1.1×5.0 mm

Rel. SNR: 1.00

SIEMENS: fl3d_rd

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	On
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	50 %
Position	L10.0 A33.9 H3.6
Orientation	Sagittal
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	26
FoV read	280 mm
FoV phase	32.8 %
Slice thickness	5.00 mm
TR	23 ms
TE	8 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	10 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	75 %
Slice resolution	54 %
Phase partial Fourier	6/8
Slice partial Fourier	6/8
Interpolation	Off
Matrix Coil Mode	CP
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off

Geometry

Multi-slice mode

Series Interleaved

Special sat. None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode None

Inline

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	On
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	Off

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Bandwidth	300 Hz/Px
Flow comp.	Yes
RF pulse type	Normal
Gradient mode	Normal
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\PC_3D_tra_fast_mip_30

TA: 3:05 PAT: Off Voxel size: 1.1×0.9×1.1 mm Rel. SNR: 1.00 SIEMENS: fl_pc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	7.7 %
Slices per slab	52
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	1.07 mm
TR	29.90 ms
TE	5.96 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	LO1-3

Contrast

Flip angle	10 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	75 %
Slice resolution	57 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Special sat.

None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Free
Encodings	2
Velocity enc. 1	30 cm/s
Velocity enc. 2	30 cm/s
Direction 1	A >> P
Direction 2	R >> L
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Bandwidth	300 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast*
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\PC_3D_tra_fast_mip_75

TA: 3:05 PAT: Off Voxel size: 1.1×0.9×1.1 mm Rel. SNR: 1.00 SIEMENS: fl_pc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	7.7 %
Slices per slab	52
FoV read	220 mm
FoV phase	100.0 %
Slice thickness	1.07 mm
TR	29.90 ms
TE	5.96 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	LO1-3

Contrast

Flip angle	10 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	75 %
Slice resolution	57 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Special sat.

None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Free
Encodings	2
Velocity enc. 1	75 cm/s
Velocity enc. 2	75 cm/s
Direction 1	A >> P
Direction 2	R >> L
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Bandwidth	300 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast*
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\TOF_3D_multi-slab brain mra

TA: 14:14 PAT: Off Voxel size: 0.5×0.5×0.5 mm Rel. SNR: 1.00 SIEMENS: fl_tof

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	6
Dist. factor	-20.00 %
Position	L0.0 A10.0 H0.0
Orientation	Transversal
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	20.0 %
Slices per slab	40
FoV read	200 mm
FoV phase	90.6 %
Slice thickness	0.50 mm
TR	20 ms
TE	3.59 ms
Averages	1
Concatenations	6
Filter	Normalize
Coil elements	LO1-3

Contrast

TD	0.000 ms
MTC	Off
Flip angle	18 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1

Resolution

Base resolution	384
Phase resolution	95 %
Slice resolution	50 %
Phase partial Fourier	Off
Slice partial Fourier	7/8
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	On
Intensity	Medium
Cut off	20
Width	4

Unfiltered images	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry

Multi-slice mode	Sequential
Series	Descending
Special sat.	Tracking H
Gap	10 mm
Thickness	40 mm

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	Head > Brain
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	L0.0 A10.0 H0.0
Orientation	Transversal
Rotation	90.00 deg
A >> P	200 mm
R >> L	182 mm
F >> H	100 mm

Physio

1st Signal/Mode	None
Dark blood	Off

Angio

TONE ramp	70 %
Flow direction	F >> H
3D centric reordering	On
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	On
MIP-Cor	On
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
--------------	----

SIEMENS MAGNETOM TrioTim syngo MR B17

Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	165 Hz/Px
Flow comp.	Yes
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\TOF_2D_obl Sagittal sinus
TA: 3:16 PAT: 2 Voxel size: 1.0×1.0×2.5 mm Rel. SNR: 1.00 SIEMENS: fl_tof

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	70
Dist. factor	-33.00 %
Position	R5.6 A27.2 H0.2
Orientation	S > C10.0
Phase enc. dir.	A >> P
Rotation	0.00 deg
Phase oversampling	0 %
FoV read	250 mm
FoV phase	100.0 %
Slice thickness	2.5 mm
TR	19 ms
TE	4.83 ms
Averages	1
Concatenations	70
Filter	Prescan Normalize, Elliptical filter
Coil elements	HEA;HEP

Contrast

TD	0.000 ms
MTC	Off
Flip angle	60 deg
Fat suppr.	None
Water suppr.	None
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Each measurement

Resolution

Base resolution	256
Phase resolution	100 %
Phase partial Fourier	Off
Interpolation	Off
PAT mode	GRAPPA
Accel. factor PE	2
Ref. lines PE	24
Matrix Coil Mode	Auto (Triple)
Reference scan mode	Integrated
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	On
Normalize	Off
B1 filter	Off
Raw filter	Off

Elliptical filter	On
Mode	Inplane
POCS	Off

Geometry

Multi-slice mode	Sequential
Series	Ascending
Sat. region 1	
Thickness	80 mm
Position	L0.0 A12.7 F145.3
Orientation	Transversal
Special sat.	None

System

Body	Off
HEP	On
HEA	On
Positioning mode	REF
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Save uncombined	Off
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R5.6 A27.2 H0.2
Orientation	S > C10.0
Rotation	0.00 deg
F >> H	250 mm
A >> P	250 mm
R >> L	119 mm

Physio

1st Signal/Mode	None
Dark blood	Off

Angio

Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	On
MIP-Cor	On
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	2D
Asymmetric echo	Allowed
Contrasts	1
Bandwidth	201 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	Yes
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\PC_3D_75_75_75

TA: 4:26 PAT: Off Voxel size: 0.9×0.9×0.3 mm Rel. SNR: 1.00 SIEMENS: fl_pc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	144
FoV read	120 mm
FoV phase	75.0 %
Slice thickness	0.25 mm
TR	38.30 ms
TE	5.57 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	LO1-3

Contrast

Flip angle	15 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	50 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Special sat.

None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default

Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Single vel.
Encodings	3
Velocity enc.	75 cm/s
Direction 1	Through plane
Direction 2	R >> L
Direction 3	A >> P
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Bandwidth	300 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast*
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\PC_3D_30_30_30_fromprev
TA: 4:46 PAT: Off Voxel size: 0.9×0.9×0.3 mm Rel. SNR: 1.00 SIEMENS: fl_pc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	144
FoV read	120 mm
FoV phase	75.0 %
Slice thickness	0.25 mm
TR	41.20 ms
TE	6.29 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	LO1-3

Contrast

Flip angle	15 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Slice resolution	50 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Special sat.

None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Single vel.
Encodings	3
Velocity enc.	30 cm/s
Direction 1	Through plane
Direction 2	R >> L
Direction 3	A >> P
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	On
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Bandwidth	300 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast*
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\PC_3D_sag_fast_mip

TA: 1:54 PAT: Off Voxel size: 1.0×0.8×3.0 mm Rel. SNR: 1.00 SIEMENS: fl_pc

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	On
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slab group 1	
Slabs	1
Dist. factor	20 %
Position	R16.9 A48.6 H11.2
Orientation	T > C-37.1
Phase enc. dir.	R >> L
Rotation	90.00 deg
Phase oversampling	0 %
Slice oversampling	6.7 %
Slices per slab	30
FoV read	200 mm
FoV phase	100.0 %
Slice thickness	3.00 mm
TR	32.45 ms
TE	6.77 ms
Averages	1
Concatenations	1
Filter	Elliptical filter
Coil elements	LO1-3

Contrast

Flip angle	10 deg
Averaging mode	Short term
Reconstruction	Magnitude
Measurements	1
Multiple series	Off

Resolution

Base resolution	256
Phase resolution	75 %
Slice resolution	56 %
Phase partial Fourier	Off
Interpolation	On
PAT mode	None
Matrix Coil Mode	Auto (CP)
Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off
Raw filter	Off
Elliptical filter	On
Mode	Inplane

Geometry

Multi-slice mode	Sequential
Series	Interleaved

Special sat.

None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
AutoAlign	---
Auto Coil Select	Default
Shim mode	Tune up
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	Isocenter
Orientation	Transversal
Rotation	0.00 deg
R >> L	350 mm
A >> P	263 mm
F >> H	350 mm

Physio

1st Signal/Mode	None
Segments	1

Angio

Flow mode	Free
Encodings	2
Velocity enc. 1	30 cm/s
Velocity enc. 2	20 cm/s
Direction 1	A >> P
Direction 2	R >> L
Rephased images	On
Magnitude images	Off
Magnitude sum	On
Phase images	Off
Subtract	Off
Std-Dev-Sag	Off
Std-Dev-Cor	Off
Std-Dev-Tra	Off
Std-Dev-Time	Off
MIP-Sag	On
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Sequence

Introduction	On
Dimension	3D
Elliptical scanning	Off
Asymmetric echo	Weak
Contrasts	1
Bandwidth	300 Hz/Px

SIEMENS MAGNETOM TrioTim syngo MR B17

Flow comp.	No
RF pulse type	Normal
Gradient mode	Fast
RF spoiling	On

SIEMENS MAGNETOM TrioTim syngo MR B17

\\USER\McLean Research\Frederick\Lia coil test\ep2d_bold_MGH_tb

TA: 1.5 s PAT: Off Voxel size: 1.2×1.2×1.8 mm Rel. SNR: 1.00 ! USER: ep2d_bold_MGH_tb

Properties

Prio Recon	Off
Before measurement	
After measurement	
Load to viewer	On
Inline movie	Off
Auto store images	On
Load to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Start measurement without further preparation	On
Wait for user to start	Off
Start measurements	single

Routine

Slice group 1	
Slices	19
Dist. factor	0 %
Position	R16.9 A48.6 H11.9
Orientation	T > C-37.1
Phase enc. dir.	P >> A
Rotation	180.00 deg
Phase oversampling	0 %
FoV read	150 mm
FoV phase	50.0 %
Slice thickness	1.8 mm
TR	1500 ms
TE	30 ms
Averages	1
Concatenations	1
Filter	Prescan Normalize
Coil elements	LO1-3

Contrast

MTC	Off
Flip angle	60 deg
Fat suppr.	Fat sat.
Averaging mode	Long term
Reconstruction	Magnitude
Measurements	1
Delay in TR	0 ms
Multiple series	Off

Resolution

Base resolution	128
Phase resolution	100 %
Phase partial Fourier	6/8
Interpolation	Off
PAT mode	None
Matrix Coil Mode	Auto (CP)
Distortion Corr.	Off
Unfiltered images	Off
Prescan Normalize	On
Raw filter	On
Elliptical filter	Off
Hamming	Off

Geometry

Multi-slice mode	Interleaved
Series	Interleaved
Sat. region 1	

Thickness	107 mm
Position	L0.0 P32.3 F24.8
Orientation	C > T37.4
Special sat.	None

System

Body	Off
LO1	On
LO2	On
LO3	On
Positioning mode	FIX
Table position	H
Table position	0 mm
MSMA	S - C - T
Sagittal	L >> R
Coronal	P >> A
Transversal	F >> H
Coil Combine Mode	Sum of Squares
AutoAlign	---
Auto Coil Select	Default
Shim mode	Standard
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Silicone	Off
? Ref. amplitude 1H	0.000 V
Adjustment Tolerance	Auto
Adjust volume	
Position	R16.9 A48.6 H11.9
Orientation	T > C-37.1
Rotation	180.00 deg
R >> L	150 mm
A >> P	75 mm
F >> H	35 mm

Physio

1st Signal/Mode	None
-----------------	------

BOLD

GLM Statistics	Off
Dynamic t-maps	Off
Starting ignore meas	0
Ignore after transition	0
Model transition states	On
Temp. highpass filter	On
Threshold	4.00
Paradigm size	20
Meas[1]	Baseline
Meas[2]	Baseline
Meas[3]	Baseline
Meas[4]	Baseline
Meas[5]	Baseline
Meas[6]	Baseline
Meas[7]	Baseline
Meas[8]	Baseline
Meas[9]	Baseline
Meas[10]	Baseline
Meas[11]	Active
Meas[12]	Active
Meas[13]	Active
Meas[14]	Active
Meas[15]	Active
Meas[16]	Active
Meas[17]	Active
Meas[18]	Active

SIEMENS MAGNETOM TrioTim syngo MR B17

Meas[19]	Active
Meas[20]	Active
Motion correction	On
Interpolation	3D-K-space
Spatial filter	Off

Sequence

Introduction	Off
Bandwidth	1260 Hz/Px
Free echo spacing	Off
Echo spacing	0.92 ms
EPI factor	64
RF pulse type	Normal
Gradient mode	Fast

Table of contents

\\USER

McLean Research
Frederick

Lia coil test

Coil is Lia Olympic
 Head First Supine
 localizer
 Make sure this scan is done with the BC
 AAHScout
 me_mprage
 epi_slicepositioning
 ep2d_bold_highres
 Five adjacent stacks - start from bottom and go up
 ep2d_bold_highres_shorttr_stack
 ep2d_bold_highres_shorttr_stack
 ep2d_bold_highres_shorttr_stack
 ep2d_bold_highres_shorttr_stack
 ep2d_bold_highres_shorttr_stack
 gre_field_mapping
 Run for Phase Contrast vessels
 vessels_head
 PC_3D_tra_fast_mip_30
 PC_3D_tra_fast_mip_75
 TOFs
 TOF_3D_multi-slab brain mra
 DO NOT SCAN PAST HERE
 TOF_2D_obl Sagittal sinus
 PC_3D_75_75_75
 PC_3D_30_30_30_fromprev
 PC_3D_sag_fast_mip
 ep2d_bold_MGH_tb